


# Low Level RF System of the LIGHT Proton Therapy Linac

Dario Soriano

AVO-ADAM RF department

12/05/2022



## Democratising Proton Therapy

*The following presentation of the AVO's LIGHT® Proton Therapy Solution is part of our Development roadmap and is subject to conformity assessment(s) by AVO's Notified Body as well as 510(k) clearance by the USA-FDA*

FORM-01180-C

# Disclaimer

This presentation may contain certain projections and other forward-looking statements with respect to the financial condition, results of operations, businesses and prospects of **Advanced Oncotherapy plc** (“**AVO**”, “**Advanced Oncotherapy**“ or the “**Company**”). These statements are based on current expectations and involve risk and uncertainty because they relate to events and depend upon circumstances that may or may not occur in the future. There are a number of factors which could cause actual results or developments to differ materially from those expressed or implied by these forward-looking statements. Any of the assumptions underlying these forward-looking statements could prove inaccurate or incorrect and therefore any results contemplated in the forward-looking statements may not actually be achieved. Nothing contained in this presentation should be construed as a profit forecast or profit estimate. Investors or other recipients are cautioned not to place undue reliance on any forward-looking statements contained herein. Advanced Oncotherapy undertakes no obligation to update or revise (publicly or otherwise) any forward-looking statement, whether as a result of new information, future events or other circumstances. Neither this presentation nor any verbal communication shall constitute an invitation or inducement to any person to subscribe for or otherwise acquire securities in Advanced Oncotherapy.

# Agenda

- [AVO-ADAM](#)
- [LIGHT - RF system](#)
- [LIGHT - LLRF system description](#)
- [LIGHT - LLRF features](#)
- [Highlights of Libera LLRF](#)

# AVO-ADAM

**2007** - Foundation of ADAM (Application of Detectors and Accelerators to Medicine), as a spin-off from CERN (the European Organisation for Nuclear Research).

**2013** - Advanced Oncotherapy's acquire ADAM.

**2018** – First LIGHT prototype in Geneva able to accelerate 52 MeV protons with energy modulation at 200 Hz.

**2022** – LIGHT Beam Production System: 230 MeV beam commissioning at Daresbury Integration V&V Site (DIS) ongoing

The LIGHT linear accelerator design by AVO-ADAM offers a modular compact solution for precise control of the treatment dose delivery, both position and energy wise. Proton energy can be modulated at up to 200 Hz in a range from 70 to 230 MeV by varying the gradient of the accelerating structures [\*]

[\*] B.B. Baricevic, A. Bardorfer, R. Cerne, G. De Michele, and Ye. Ivanisenko, “Light Proton Therapy Linac LLRF System Development”, in *Proc. LINAC'18*, Beijing, China, Sep. 2018, pp. 171-173. doi:10.18429/JACoW-LINAC2018-MOP0081



# LIGHT - RF system

Recent advances at DIS, click on the link or scan QR:  
<https://www.youtube.com/watch?v=07TFUL5DzMU>




## LIGHT accelerating cavities:

- 1 RFQ (750 MHz)
- 4 SCDLs (3 GHz)
- 15 CCLs (3 GHz)

## RF generation & amplification

- 14 High Power Stations, providing 200 Hz, 5  $\mu$ s, RF pulses up to 7.5 MW
- 14 LLRF boxes - Libera LLRF by Instrumentation Technologies



# LIGHT - LLRF system description

## 2 types of LLRF boxes:

- a. 749.48 MHz
  - Digital Processor + Splitter Unit
- b. 2997.92 MHz with 1 or 2 ADC boards
  - Front End: temperature stabilization of the analog PCBs inside the LLRF analog front-end unit
  - Digital Processor: where the IF signals are processed by the FPGAs and a drive signal is generated



## Other LLRF modules

- Reference Master Oscillator + Distribution amplifier
  - CW reference at 2997.92 MHz to distribute the signal to each LLRF unit.
- Interlock unit
  - Distributes the interlock signals from the accelerator control system to each LLRF box.
- Trigger Synchronization unit



One of the four LLRF racks in DIS during pre-installation SAT period


# LIGHT - LLRF features (1/2)

- ✓ **Feedback control** - Amplitude and phase regulation
- ✓ **RF output stability** - Amplitude and phase
  - Requirement: better than 0.15 % in amplitude and 0.15 deg in phase
  - Tested, values ~ 0.01 % and 0.01 deg RMS at full range.
- ✓ **Decay Analysis**
  - RF signal analyzed to provide a virtual frequency correction (frequency is fixed for beam acceleration).
  - This can be translated into a temperature correction to keep resonance of the cavity by changing water temperature.
  - Via Modbus RTU communication protocol
  - Tested standalone, ongoing work to incorporate into the full LINAC.
- ✓ **Pulse shapes**
  - Possibility to define second order polynomials to shape the output pulse.
  - Flatten RF pulse after HPS amplification.



RFQ voltage (purple) during a test where cooling (white) was oscillating excessively. Once feedback control loop (cyan) was activated, stability improved.

Shape of the RF pulse in the RF network. Uncompensated.



# LIGHT - LLRF features (2/2)

## ✓ Breakdown detection

- Detects when reflected RF power rises over a threshold.
- Counts number of events per unit time.
- Stops RF output if events are above allowed limit.
- Cavity protection from arcing discharges.

## ✓ Front End Controller

- Data Streaming Interface
  - Uni-directional interface to receive pulse related measurements..
- Slow Control Interface
  - Bi-directional interface to provide control, configuration and calibration settings to the system and to collect acquisition data that has update intervals in the 100 milliseconds.
- Real Time Interface
  - Bi-directional deterministic request-response interface to configure pulse-related settings.
  - All messages and responses within 5 ms.
- Trigger Interface
  - Uni-directional interface to send trigger signals to the LLRF Unit.



# Highlights of Libera LLRF

## - **Reliability**

- Barely zero downtime related to LLRF units during beam commissioning in DIS

## - **Configurability**

- Many parameters and options to adjust the LLRF output / control. Possible to tailor it to the application needed.
- Possibility to create routines easily and ease acceptance testing or calibration.

## - **Graphical User Interface**

- For testing purposes. It was quite handy before our control system was developed, allows to quickly see what you are doing and have the knobs one click away.

## - **Support from I-tech**

- Regular meetings > For updates on the project and purchases
- Technical meetings > For specific questions, issues or ideas.
- Email exchanges > Detailed descriptions supported by the LLRF documentation (>100 pages)

AVO-ADAM and I-tech collaboration continues for:

- Further improvement of the LLRF firmware
  - Production spare-extra LLRF units

### **GIOVANNI DE MICHELE**

HEAD OF RF GROUP & DEPUTY TECH. DIRECTOR (AVO-ADAM)

Instrumentation Technologies is a key supplier and an industrial partner providing the LLRF controls and BPM electronics for the LIGHT system. Via the Libera platforms, with ad-hoc features design, I-Tech has met the requirements of the LIGHT system. We are pleased to collaborate with such a professional, open and flexible team that has helped us enhance our customer experience.

**Thank you very much !**

**Happy to take questions**